

PRODUCT INFORMATION

pUC19 DNA/Mspl (HpaII) Marker, 23

#SM0221 **50 µg**
(for 100 applications)

Lot: —

Concentration: 0.5 µg/µL
Supplied with: 1 mL 6X DNA Loading Dye

Store at -20°C

In total 2 vials.

www.thermoscientific.com/onebio

Description

pUC19 DNA was completely digested with Mspl, purified and dissolved in a storage buffer.

The DNA Marker contains the following 13 discrete fragments (in base pairs): 501, 489, 404, 331, 242, 190, 147, 111, 110, 67, 34, 34, 26.

Storage Buffer

10 mM Tris-HCl (pH 7.6), 1 mM EDTA.

6X DNA Loading Dye

10 mM Tris-HCl (pH 7.6), 0.03% bromophenol blue, 0.03% xylene cyanol FF, 60% glycerol and 60 mM EDTA.

CERTIFICATE OF ANALYSIS

Well-defined bands are formed during agarose gel electrophoresis.

The DNA concentration is determined spectrophotometrically.

The absence of nucleases is confirmed by a direct nuclease activity assay.

Quality authorized by:



Jurgita Zilinskiene


I. Loading on agarose gel:

- prepare the marker before loading:
 - 1 μL (0.5 μg) of the marker,
 - 1 μL of 6X DNA Loading Dye,
 - 4 μL of deionized water;
- vortex gently just prior to use;
- do not heat before loading;
- apply the prepared amount (6 μL) of the marker on a 5 mm lane of agarose gel;
- following electrophoretic separation on gel, visualize the DNA bands by ethidium bromide staining.

Important Notes

- One vial (50 μg) is sufficient for \sim 100 applications.
- Use 0.1 μg (0.2 μL) of the marker (before dilution) per 1 mm of an agarose gel lane width.
- For DNA band visualization with SYBR[®] Green, GelRed and other intercalating dyes, do not add the dyes into the sample, use gel staining after electrophoresis or include dyes into agarose gel to avoid aberrant DNA migration.

pUC19 DNA/MspI (HpaII) Marker, 23



*The 501 and 489 bp bands migrate anomalously (1, 2, 3).
26 bp fragment is not visible and comprises 1.0%

(continued on back page)

II. Loading on polyacrylamide gel (1, 2, 3):

- prepare the marker before loading:
 - 2 μL (1 μg) of the marker,
 - 0.5 μL of 6X DNA Loading Dye,
 - 0.5 μL of deionized water;
- vortex gently just prior to use;
- do not heat before loading;
- apply the prepared amount (3 μL) of the marker on a 5 mm lane of polyacrylamide gel;
- following electrophoretic separation on gel, visualize the DNA bands by ethidium bromide staining.

Note

- One vial (50 μg) is sufficient for ~50 applications.
- Use 0.2 μg (0.4 μL) of the marker per 1 mm of a polyacrylamide gel lane width.

References

1. Stellwagen, N.C., Anomalous electrophoresis of deoxyribonucleic acid restriction fragments on polyacrylamide gels, *Biochemistry*, 22, 6186-6193, 1983.
2. Lane, D., et al., Use of gel retardation to analyze protein – nucleic acid interactions, *Microbiological Reviews*, 56, 509-528, 1992.
3. Stellwagen, N.C., Conformational isomers of curved DNA molecules can be observed by polyacrylamide gel electrophoresis, *Electrophoresis*, 21, 2327-2334, 2000.

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.thermoscientific.com/onebio for Material Safety Data Sheet of the product.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. SYBR is a registered trademark of Molecular Probes, Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.